



**THE  
BINDING  
SITE**

# **FREELITE™**

*For Use On*  
**Bayer ADVIA® 1650**

## **Highly Specific Reagents for the Measurement of Immunoglobulin Free Light Chains in Serum**

### **SENSITIVE**

Latex enhanced reagents give previously unattainable levels of sensitivity.  
<1.0mg/L for serum **free kappa** and **free lambda**<sup>1</sup>.

### **SPECIFIC**

Specificity confirmed against whole immunoglobulin preparations. Minimal cross reactivity allows true free light chain levels to be precisely measured in serum despite myelomatous levels of whole immunoglobulin<sup>2</sup>.

### **SIMPLE**

Full application data and detailed instructions on programming are available. All reagents necessary for running the assays, additional to standard system buffers, are supplied in the kits.

### **RELIABLE**

Reference ranges have been established and published in *Clinical Chemistry*<sup>3</sup>. Many of the world's leading reference laboratories are now routinely offering Freelite assays as part of their Myeloma investigation panel.



#### **References:**

1. Serum Free Light Chain Analysis. AR Bradwell, GP Mead, HD Carr-Smith 2005.
2. Highly Sensitive, Automated Immunoassay for Immunoglobulin Free Light Chains in Serum and Urine. Arthur R. Bradwell, Hugh D. Carr-Smith, Graham P. Mead, Lian X. Tang, Paul J. Showell, Mark T. Drayson and Roger Drew. *Clin. Chem*, 2001;47:4;673-680.
3. Serum Reference Intervals and Diagnostic Ranges for Free  $\kappa$  and Free  $\lambda$  Immunoglobulin Light Chains: Relative Sensitivity for Detection of Monoclonal Light Chains. Jerry A. Katzmann, Raynell J. Clark, Roshini S. Abraham, Sandra Bryant, James F. Lymp, Arthur R. Bradwell and Robert A. Kyle. *Clin. Chem*, 2002;48:9;1437-1444.

Freelite™ is a trademark of The Binding Site Ltd, Birmingham, UK.  
Bayer ADVIA® is a registered trademark of Bayer Corporation, New York, USA.

# FREELITE™ for use on Bayer ADVIA® 1650

	KAPPA	LAMBDA
<b>NORMAL ADULT SERUM RANGE</b> (95 percentiles)	3.3 - 19.4 mg/L	5.7 - 26.3 mg/L
<b>MEASURING RANGE</b> Standard sample dilution	3.74 - 56.2 mg/L	7.0 - 93.4 mg/L
<b>SENSITIVITY</b> Minimum serum sample dilution (1/1)	0.75 mg/L	0.88 mg/L
<b>STANDARD SERUM SAMPLE DILUTION</b>	1/5	1/8
<b>MINIMUM SERUM SAMPLE DILUTION</b>	1/1	1/1
Sample volume per test (diluted)	5 microlitres	3 microlitres
Minimum sample volume	50 microlitres	50 microlitres
Reagent volume per test	55 microlitres	55 microlitres
Dead volume of reagent per vial	2000 microlitres	2000 microlitres
<b>INTRA ASSAY PRECISION</b>		
Level 1	6.51%	6.28%
Level 2	2.72%	6.42%
Level 3	2.34%	1.81%
<b>INTER ASSAY PRECISION</b>		
Level 1	4.98%	7.07%
Level 2	1.91%	2.90%
Level 3	2.59%	2.96%
<b>OPEN VIAL STABILITY</b>		
Latex reagent - Liquid ready to use	3 months	3 months
Calibrators & Controls - Liquid ready to use	3 months	3 months
<b>KIT CONTENTS</b>		
Latex reagent	2 x 4.0 mL	2 x 4.0 mL
Calibrator	2 x 1.5 mL	2 x 1.5 mL
Low Control	1 x 1.0 mL	1 x 1.0 mL
High Control	1 x 1.0 mL	1 x 1.0 mL
Supplementary Reagent	1 x 25 mL	1 x 25 mL

ANALYSER	DESCRIPTION	PACK	CODE
Bayer ADVIA® 1650	Freelite Kappa Kit	2 x 50 test	LK016.B
	Freelite Lambda Kit	2 x 50 test	LK018.B

All kits are FDA approved for *in vitro* diagnostic use. See product insert for further information. Kits are CE marked for most European countries, please contact us for the latest information.

**MKG286.E**

Nov 05

THE BINDING SITE LIMITED

P.O. Box 11712, Birmingham, B14 4ZB, UK.

Tel: +44 (0)121 436 1000 Fax: +44 (0)121 430 7061 info@bindingsite.co.uk

www.freelite.co.uk www.bindingsite.co.uk

