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Auteur	N	Détail
van Rhee <i>et al</i> ²¹	301	"High baseline SFLC levels were a reflection of higher tumor burden, higher degree of disease aggressiveness and light-chain-only MM with its greater propensity for renal failure."
Kyrtsolis <i>et al</i> ²²	94	"The 5-year disease-specific survival was 82% and 30% in patients with [serum free light chain ratio] sFLCR lower than and equal or greater than the median, respectively (P=0.0001). sFLCR was an independent prognostic factor."
Snozek <i>et al</i> ²³	790	"These findings suggest that the serum FLC ratio at initial diagnosis is an important predictor of prognosis in myeloma, and can be incorporated into the ISS for improved risk stratification."
Kyrtsolis <i>et al</i> ²⁴	214	"In conclusion, baseline sFLCR appears to be an easily determined powerful, independent and very promising novel prognostic factor for survival in patients with newly diagnosed MM."
Dispenzieri <i>et al</i> ²⁵	273	"An abnormal ratio proved to be an independent predictor of adverse outcome - in the case of SMM - progression to active MM."

"Unlike baseline and follow-up analyses of serum and urine M-proteins, high sFLC levels at baseline-reflecting more aggressive disease- and steeper reductions after therapy identified patients with inferior survival"

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Références

1. Bakshi *et al.* Serum Free Light Chain (FLC) Measurement Can Aid Capillary Zone Electrophoresis in Detecting Subtle FLC-Producing M Proteins. *Am J Clin Pathol* 2005; **124**: 214-218
2. Abadie *et al.* Assessment of Serum Free Light Chain Assays for Plasma Cell Disorder Screening in a Veterans Affairs Population. *Ann Clin & Lab Sci* 2006; **36**:157-162
3. Hill *et al.* Serum Free Light Chains: An Alternative to the Urine Bence Jones Proteins Screening Test for Monoclonal Gammopathies. *Clin Chem* 2006; **52**:1743-1748
4. Smith *et al.* Guidelines on the diagnosis and management of multiple myeloma. *Br J Haematol* 2005; **132**:410-451
5. Katzmam *et al.* Elimination of the Need for Urine Studies in the Screening Algorithm for Monoclonal Gammopathies by Using Serum Immunofixation and Free Light Chain Assays. *Mayo Clin Proc* 2006; **81**:1575-1578
6. Mead *et al.* Serum free light chains for monitoring multiple myeloma *Br J Haematol* 2004; **126**:348-354
7. Bradwell *et al.* Serum test for assessment of patients with Bence Jones myeloma. *Lancet* 2003; **361**:489-491
8. Tate *et al.* Analytical performance of serum free light-chain assay during monitoring of patients with monoclonal light-chain diseases. *Clin Chim Acta* 2007; **376**:30-36
9. Abraham *et al.* Correlation of Serum Immunoglobulin Free Light Chain Quantification with Urinary Bence Jones Protein in Light Chain Myeloma. *Clin Chem* 2002; **48**:655-657
10. Katzmam *et al.* Diagnostic Performance of Quantitative κ and λ Free Light Chain Assays in Clinical Practice. *Clin Chem* 2005; **51**:878-881
11. Drayson *et al.* Serum free light-chain measurements for identifying and monitoring patients with nonsecretory multiple myeloma. *Blood* 2001; **97**: 2900-2902
12. Lachmann *et al.* Outcome in systemic AL amyloidosis in relation to changes in concentration of circulating free immunoglobulin light chains following chemotherapy. *Br J Haematol* 2003; **122**:78-84
13. Abraham *et al.* Quantitative Analysis of Serum Free Light Chains, A New Marker for the Diagnostic Evaluation of Primary Systemic Amyloidosis. *Am J Clin Pathol* 2003; **119**:274-278
14. Katzmam *et al.* Serum Reference Intervals and Diagnostic Ranges for Free κ and Free λ Immunoglobulin Light Chains: Relative Sensitivity for Detection of Monoclonal Light Chains. *Clin Chem* 2002; **48**:1437-1444
15. Dingli *et al.* Immunoglobulin free light chains and solitary plasmacytoma of bone. *Blood* 2006; **108**:1979-1983
16. Rajkumar *et al.* Serum free light chain ratio is an independent risk factor for progression in monoclonal gammopathy of undetermined significance. *Blood* 2005; **106**:812-817
17. Hajek *et al.* Free Light Chain Assays for Early Detection of Resistance to Bortezomib-Based Regimens. *Haematologica* 2007;**92**:257a
18. Mösbauer *et al.* Monitoring serum free light chains in patients with multiple myeloma who achieved negative immunofixation after allogenic stem cell transplantation. *Haematologica* 2007; **92**:275-276
19. Dawson *et al.* Extramedullary relapse of multiple myeloma associated with a shift in secretion from intact immunoglobulin to light chains. *Haematologica* 2007; **92**:143-144
20. Mead *et al.* Incidence of light chain escape in myeloma patients at relapse. *Br J Haematol* 2008; **141**:98a
21. van Rhee *et al.* High serum-free light chain levels and their rapid reduction in response to therapy define an aggressive multiple myeloma subtype with poor prognosis. *Blood* 2007; **110**:827-832
22. Kyrtonis *et al.* Prognostic value of serum free light chain ratio at diagnosis in multiple myeloma. *Br J Haematol* 2007; **137**:240-243
23. Snozek *et al.* Prognostic value of the serum free light chain ratio in newly diagnosed myeloma: proposed incorporation into the international staging system. *Leukemia* advance online publication, 3 July 2008;doi:10.1038/leu.2008.171
24. Kyrtonis *et al.* The Addition of sFLCR Improves ISS Prognostication in Multiple Myeloma (MM). *Blood* 2007; **110**:1490a
25. Dispenzieri *et al.* Immunoglobulin free light chain ratio is an independent risk factor for progression of smoldering (asymptomatic) multiple myeloma. *Blood* 2008; **111**:785-789
26. Gertz *et al.* Definition of Organ Involvement and Treatment Response in Immunoglobulin Light Chain Amyloidosis (AL): A Consensus Opinion From The 10th International Symposium on Amyloid and Amyloidosis. *Am J Hematol* 2005; **79**:319-328
27. Bird *et al.* Guidelines on the diagnosis and management of AL amyloidosis. *Br J Haematol* 2004; **125**:681-700
28. National Academy of Clinical Biochemistry. Laboratory Medicine Practice Guidelines, Practical Guidelines And Recommendations For Use Of Tumor Markers In The Clinic. www.aacc.org/members/nacb/LMPG
29. Durie *et al.* International uniform response criteria for multiple myeloma. *Leukemia* 2006; **20**:1467-1473

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Cette recherche a été initialement publiée dans *Blood*: Rajkumar *et al*. Serum free light chain ratio is an independent risk factor for progression in monoclonal gammopathy of undetermined significance. *Blood* 2005; **106**:812-817 ©2005 The American Society of Hematology.
Germany & Austria: BINDING SITE GmbH